

DeVilbiss® Atomizers and Tips



Use convenient, single-use RhinoGuard™ Tip Covers on your DeVilbiss products to reduce cross contamination, increase comfort for your patients and eliminate the need to decontaminate products over and over.



ATOMIZER

Ultra-fine particle atomization of most medical solutions

Straight flow tip is easily disassembled for cleaning

Suitable for applying fine, medicated spray to nasal and oral cavities, or oral vaccines

May be used with compressed air, the No.1 bulb or No. 630 cut-off assembly

Model DV-286

Atomizer, clear glass with metal top
1 per package/non-sterile



ATOMIZER

Designed for quick change and replacement of spray tube

Tight fitting "O" ring ensures consistent spray

Metal spray tube may be sterilized after each use

Adjustable tip permits up or down spray for aqueous solutions of heavier viscosity

Operated by the No.630 cutoff assembly or a No. 1 bulb

Model DV-163

Atomizer, amber glass with metal top
1 per package/non-sterile

RHINO-GUARD ATOMIZER TIP PROTECT

SP-30010 100 per box/non-sterile

SP-30250 250 per box/non-sterile



5860-C Hollister Avenue
Santa Barbara, CA 93117 U.S.A.

Telephone: 805 967-1153
FAX: 805 967-1295

Website: www.tecfen.com
email: tecfen@tecfen.com

