

Eavey E3 Ear Wrap



SB-41400 Single Ear

Surgeon Features

- Time and money savings
- Adjustable
- One size fits all
- Use on either ear
- Change quick and easily
- Easily repositioned
- Easily replaced
- Used in combination with familiar cotton and gauze absorbent materials

Patient Features

- Extremely comfortable
- Safe sleeping
- Easily repositioned
- Easily replaced at home following surgical instructions

Surgical External Post-Operative Dressing

**Developed in conjunction with Roland D. Eavey, M.D.,
Massachusetts Eye & Ear Infirmary**

Instructions:

In the operating room, apply the conventional absorbent material such as cotton and gauze to the operated ear. Place the widest part over the operated ear and absorbent materials. Wrap the band around the head and use the velcro to secure. The wrap is placed at the same level one would have placed a conventional dressing.



5860-C Hollister Avenue
Santa Barbara, CA 93117 U.S.A



Tel: (805) 967 1153
Fax: (805) 967 1295



Home: www.tecfe.com
Email: tecfe@tecfe.com