

Trephine Blades, Traditional, Non-Vacuum, and Vacuum



Traditional Corneal Trephines

- Long and short styles
- Compatible with most punch systems
- Reproducible results
- Range 3.00 to 18.00 mm long and 6.00 to 9.00 mm short
- Eye Bank sizes available

LONG TREPHINES

TE-330300, 3.00mm Diameter	TE-330900, 9.00mm Diameter
TE-330400, 4.00mm Diameter	TE-330950, 9.50mm Diameter
TE-330450, 4.50mm Diameter	TE-331000, 10.00mm Diameter
TE-330500, 5.00mm Diameter	TE-331050, 10.50mm Diameter
TE-330550, 5.50mm Diameter	TE-331100, 11.00mm Diameter
TE-330575, 5.75mm Diameter	TE-331150, 11.50mm Diameter
TE-330600, 6.00mm Diameter	TE-331200, 12.00mm Diameter
TE-330625, 6.25mm Diameter	TE-331250, 12.50mm Diameter
TE-330650, 6.50mm Diameter	TE-331300, 13.00mm Diameter
TE-330675, 6.75mm Diameter	TE-331350, 13.50mm Diameter
TE-330700, 7.00mm Diameter	TE-331400, 14.00mm Diameter
TE-330725, 7.25mm Diameter	TE-331450, 14.50mm Diameter
TE-330750, 7.50mm Diameter	TE-331500, 15.00mm Diameter
TE-330775, 7.75mm Diameter	TE-331550, 15.50mm Diameter
TE-330785, 7.85mm Diameter	TE-331600, 16.00mm Diameter
TE-330800, 8.00mm Diameter	TE-331650, 16.50mm Diameter
TE-330825, 8.25mm Diameter	TE-331700, 17.00mm Diameter
TE-330850, 8.50mm Diameter	TE-331800, 18.00mm Diameter
TE-330875, 8.75mm Diameter	Packaged sterile, one per package

SHORT TREPHINES

TE-310600 , 6.00mm Diameter
TE-310625 , 6.25mm Diameter
TE-310650 , 6.50mm Diameter
TE-310675 , 6.75mm Diameter
TE-310700 , 7.00mm Diameter
TE-310725 , 7.25mm Diameter
TE-310750 , 7.50mm Diameter
TE-310775 , 7.75mm Diameter
TE-310800 , 8.00mm Diameter
TE-310825 , 8.25mm Diameter
TE-310850 , 8.50mm Diameter
TE-310875 , 8.75mm Diameter
TE-310900 , 9.00mm Diameter

Packaged sterile, one per package

TREPHINE System & Components



TE-330010
Corneal Punch System

- 7-9.0mm Punch Piston
- Piston Holder
- Base Plate



Punch Piston

TE-330005 5-7.0mm Punch Piston
TE-330007 7-9.0mm Punch Piston
TE-330009 9-11.0mm Punch Piston
TE-330001 Disposable Corneal Block (sterile)



5385 Hollister Avenue, No. 216
Santa Barbara, CA 93111 U.S.A.



Telephone: 805 967-1153
FAX: 805 967-1295



Website: www.tecfenmedical.com
email: tecfen@tecfen.com

NON-VACUUM TREPHINE PUNCHES

TE-340600, 6.00mm Diameter
TE-340625, 6.25mm Diameter
TE-340650, 6.50mm Diameter
TE-340675, 6.75mm Diameter
TE-340700, 7.00mm Diameter
TE-340725, 7.25mm Diameter
TE-340750, 7.50mm Diameter
TE-340775, 7.75mm Diameter
TE-340800, 8.00mm Diameter
TE-330825, 8.25mm Diameter
TE-330850, 8.50mm Diameter
TE-330875, 8.75mm Diameter
TE-330900, 9.00mm Diameter
TE-330950, 9.50mm Diameter

- Posts ensure proper alignment
- Reliable performance
- Complete system Range 6.00 to 9.50mm



Packaged sterile, one per package

VACUUM TREPHINES

- Unsurpassed sharpness and consistency
- Access to all eyes, regardless of depth
- Non-slip rings on finger grips assist with handling when gloves wet
- Range 6.00 to 9.50 mm

TE-350600, 6.00mm Diameter
TE-350625, 6.25mm Diameter
TE-350650, 6.50mm Diameter
TE-350675, 6.75mm Diameter
TE-350700, 7.00mm Diameter
TE-350725, 7.25mm Diameter
TE-350750, 7.50mm Diameter
TE-350775, 7.75mm Diameter
TE-350800, 8.00mm Diameter
TE-350825, 8.25mm Diameter
TE-350850, 8.50mm Diameter
TE-350875, 8.75mm Diameter
TE-350900, 9.00mm Diameter
TE-350950, 9.50mm Diameter

Packaged sterile, one per package



Non-slip finger grips

